



4320 Forest Park Ave
Suite 303
Saint Louis, MO | 63108
+1 (314) 833-9764

E. COLI HOST CELL PROTEIN ELISA KIT

Complete kit for the determination of E. coli host cell protein contamination in bulk products expressed in E. coli expression systems.

KIT INCLUDES

Coated 96-Well Strip Plate
E. coli Protein Standard
5x Dilution Buffer
10x PBS-T
Reporting Antibody
Streptavidin-HRP Conjugate
TMB Substrate
Stop Solution
Plate Sealer

PRODUCT MANUAL HCP-002

E. COLI HOST CELL PROTEIN ELISA KIT

ASSAY PRINCIPLE

The E. coli Host Cell Protein (HCP) ELISA kit is designed to quantitatively measure host cell protein contamination in bulk products expressed in E. coli expression systems. Please read the complete kit insert before performing this assay. A series of E. coli protein standards are made from a provided stock solution and used to generate a standard curve for the assay, and all unknown sample concentrations should be read off this standard curve. E. coli standards or diluted unknown samples are pipetted into the provided 96-well plate which has been precoated with anti-E. coli HCP antibodies to capture E. coli proteins from biologics samples. Following an incubation to allow capture of the E. coli protein by the antibodies on the plate, a second anti-E. coli HCP antibody, conjugated with biotin, is added and incubated to allow it to bind to the captured E. coli proteins. After a 45 minute incubation, the plate is washed and a Streptavidin-HRP (Horse Radish Peroxidase) conjugate is added and incubated for 30 minutes. The Streptavidin-HRP conjugate will be captured by any biotin-labeled antibody bound to the plate. Following a wash step to remove unbound conjugate, TMB substrate is added and is converted by the captured HRP to a colored product in proportion to the amount of HCP bound to the plate. After a ten minute incubation to allow color development, the reaction is stopped and the intensity of the generated color is detected in a spectrophotometer plate reader capable of measuring 450 nm wavelength. A standard curve will be generated from the E. coli protein standards and used to calculate the concentration of E. coli proteins in the unknown samples, after making suitable correction for the dilution of the sample.

SUPPLIED COMPONENTS

ENTIRE KIT MUST BE STORED AT 4°C.

Clear 96-Well Strip Plate

Clear plastic strip-well microtiter plate coated with rabbit anti-E. coli HCP IgG. Can be used as individual strips.

E. coli Protein Standard (2430 ng/ml, 600 µl)

Concentrated E. coli proteins sufficient for generating a standard curve from 810 ng/ml to 1.1 ng/ml.

5x Dilution Buffer (15 ml)

1x Dilution Buffer is used for dilution of Reporting Antibody and Streptavidin-HRP conjugate. 1x Dilution Buffer is used to dilute samples if necessary. The 10 ml of 5x concentrate should be diluted to 50 ml with 40 ml of milliQ water to achieve 1x Dilution Buffer.

10x PBS-T (30 ml)

1x PBS-T is used for wash steps. The 30 ml of 10x concentrate should be diluted to 300 ml with 270 ml of milliQ water to achieve 1x PBS-T.

Reporting Antibody (1 mg/ml, 150 µl/tube)

A biotin labeled rabbit polyclonal antibody specific for E. coli cell proteins. Immediately prior to the assay, dilute the entire 75 µl into 15 ml of 1x Dilution Buffer to give a 5 µg/ml working stock.

Streptavidin-HRP Conjugate (4 µg/ml, 375 µl/tube)

A Streptavidin-Horse Radish Peroxidase conjugate in a stabilizing solution. Immediately prior to the assay, dilute the entire 375 µl into 15 ml of 1x Dilution Buffer to give a 0.1 µg/ml working stock.

TMB Substrate (15 ml)

Use directly without dilution.

Stop Solution (15 ml)

A 1M solution of sulfuric acid. CAUSTIC. Use directly without dilution. The solution should not come in contact with skin or eyes. Take appropriate precautions when handling this reagent.

Plate Sealer

OTHER MATERIALS REQUIRED

milliQ water

Single- and multi-channel micro-pipettes with disposable tips to accurately dispense volumes 5-250 µl.

1.5 ml eppendorf tubes

Reagent reservoirs for sample addition


Colorimetric 96-well microplate reader capable of reading optical density at 450 nm.



ASSAY PROTOCOL

1. Dilute the 10x PBS-T and 5x Dilution Buffer to 1x-strength with milliQ water. Check both concentrate bottles for precipitates before proceeding and if found warm slightly in a water bath to dissolve before proceeding. The 30 ml of 10x PBS-T should be diluted to 300 ml with 270 ml milliQ water. 15 ml of 5x Dilution Buffer should be diluted to 75 ml with 60 ml of milliQ water.
2. Prepare the HCP standards by numbering eight 1.5 ml tubes and adding 1ml of Dilution Buffer to each. Cap the eighth tube, this will be the blank (0 ng/ml HCP). To tube one add 500 μ l of the provided 2430 ng/ml HCP stock and mix well, this will be the 810 ng/ml standard. Then serially dilute 500 μ l of tube one across tubes two through seven to generate the remainder of the standards, with an 8th tube serving as a blank using Dilution Buffer alone. Pipette 100 μ l of each standard and the blank into the plate.
3. Pipette 100 μ l of samples into wells in the plate. If necessary, first dilute the samples in 1x Dilution Buffer.
4. Cover plate with individual plate seal and incubate **1.5 hours** at room temperature.
5. During the above incubation, dilute the 1 mg/ml Reporting Antibody to 5 μ g/ml by adding 75 μ l to 15 ml of 1x Dilution Buffer.
6. Wash plate by emptying contents and adding 250 μ l of 1x PBS-T to each well. Empty wells again and tap the plate firmly upside down on a paper towel to fully empty well. **Repeat 1x PBS-T wash step two additional times.**
7. Pipette 100 μ l of Reporting Antibody into each well. Cover plate with the plate seal and incubate plate **45 minutes** at room temperature.
8. During the above incubation, dilute the 4 μ g/ml Streptavidin-HRP conjugate to 0.1 μ g/ml by adding the entire 375 μ l to 15 ml of 1x Dilution Buffer.
9. Wash plate by emptying contents and adding 250 μ l of 1x PBS-T to each well. Empty wells again and tap the plate firmly upside down on a paper towel to fully empty well. **Repeat 1x PBS-T wash step two additional times.**
10. Pipette 100 μ l of Streptavidin-HRP conjugate into wells. Cover plate and incubate plate **30 minutes** at room temperature.
11. Wash plate by emptying contents and adding 250 μ l of 1x PBS-T to each well. Empty wells again and tap the plate firmly upside down on a paper towel to fully empty well. **Repeat 1x PBS-T wash step two additional times.**
12. Add 100 μ l of TMB substrate to each well. Monitor color development. Generally 10 minutes time will be sufficient; incubating longer may increase the background.
13. Stop reaction by adding 100 μ l of Stop Solution to each well containing TMB when the color development within standards is sufficient.
14. Read the optical density generated from each well in a plate reader capable of reading at 450 nm. A 4-parameter fit standard curve should be generated from the E. coli protein standards and used to calculate the concentration of E. coli proteins in the unknown samples, taking into account any unknown sample dilution made.

STANDARD CURVE



Lorem ipsum

HOST CELL PROTEIN DETECTION

Canopy Biosciences Host Cell Protein Detection Kits including antibodies are to be used for impurity testing and research only. Canopy Biosciences products may not be reverse engineered without prior written approval of Canopy Biosciences, LLC. **DISCLAIMER TO THE EXTENT ALLOWED BY LAW, CANOPY BIOSCIENCES, LLC WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.**

LIMITED LABEL USE LICENSE

Canopy Biosciences Host Cell Protein Detection kit and antibody Notice to Purchaser: Customer may not isolate, extract, reverse engineer, derive, copy or separately use any component of the Product (such as, for example, any antibody component) for manufacturing of additional antibody (These Products are solely for impurity testing and protein analysis as set forth in this limited use label license). Canopy Biosciences prohibits any use of antibodies as therapeutics, diagnostics, or prophylactics ; or any sale, resale, leasing, or licensing without written approval of Canopy Biosciences. For information on obtaining additional rights to Canopy Biosciences Host Cell Protein Detection products for any use not permitted herein, please contact Canopy Biosciences at info@canopybiosciences.com.

TRADEMARKS

All trademarks are the property of Canopy Biosciences and its subsidiaries unless otherwise specified